

Red Light Eye Lift Device

Combination Red Light + Microcurrent Therapy



Manual

Please read the instructions
carefully before use

Operating Instructions

Power / Mode Button

- Power On:**
Press and hold the power button for 1–3 seconds to turn on the device. The red indicator light will turn on.
- Power Off:**
Press and hold the power button for 1–3 seconds to turn off the device. The indicator light will turn off.
- Mode Selection:**
After powering on, the device will default to the last mode used.
There are three modes (Mode 1, Mode 2, Mode 3). Press the power/mode button briefly to switch between modes.

+ / – Buttons (Intensity Control)



- After startup, the device defaults to low level. Each mode has three intensity levels: low, medium, and high.
Press the + / – button to adjust the intensity level.

Charging

- Connect the Type-C charging cable to charge the device.
 - Red light: Charging
 - Green light: Fully charged

💡 **Tip:** Fully charge the device before first use or after extended periods of inactivity.

Product Overview



Product name	RLTEYE
Input voltage	5V
Rated power	3.7W
Body material	ABS

Operating Instructions



STEP 1

Clean your face using your regular cleanser, then pat dry.



STEP 2

Apply moisturizer, eye cream, or serum to the under-eye area to keep the skin hydrated.



STEP 3

Place the eye mask comfortably on your face and select your desired mode. *Tip: For first-time use, start on a low setting and gradually increase as needed.*



STEP 4

After use, remove the eye mask and store it properly.

💡 **Tip:** Ensure the inner surface of the device sits closely against your skin for best results.

Precautions

- It is normal to feel a mild tingling sensation during use. The intensity may vary depending on skin moisture—the more hydrated the skin, the stronger the sensation.
- Individual facial structure, muscle tone, and sensitivity vary from person to person. Each unit is tested for quality before leaving the factory.
- For adult use only. Always begin at the lowest intensity setting.
- Follow the instructions carefully. Do not use the device continuously on the same area for extended periods, as this may cause skin irritation or muscle discomfort.
- Not recommended for use during pregnancy or menstruation. If use is necessary, consult a healthcare professional.
- Do not use fast chargers or non-approved charging devices.
- Charge the device before the battery is fully depleted. Avoid leaving the battery completely drained for long periods.
- After use, clean any residue from the device and keep it dry. Do not allow water or moisture to enter the product.
 - This product is not waterproof.
 - Do not rinse with water.
 - If cleaning is needed, wipe with a damp cloth or tissue.